

# Product Datasheet

## Human PLA2R1 - Ready-To-Use ELISA Kit (Colorimetric) NBP3-39349

Unit Size: 1 Kit

Storage of components varies. See protocol for specific instructions.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-39349](http://www.novusbio.com/NBP3-39349)

Updated 10/16/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-39349](http://www.novusbio.com/reviews/destination/NBP3-39349)




**NBP3-39349****Human PLA2R1 - Ready-To-Use ELISA Kit (Colorimetric)**

| <b>Product Information</b> |                                                                                                          |
|----------------------------|----------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>           | 1 Kit                                                                                                    |
| <b>Concentration</b>       | Concentration is not relevant for this product. Please see the protocols for proper use of this product. |
| <b>Storage</b>             | Storage of components varies. See protocol for specific instructions.                                    |

| <b>Product Description</b>  |                                                                                                                                                                                                                 |
|-----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>          | The Ready-To-Use ELISA kit offers pre-diluted detection reagents and a shorter experimental time.<br>Assay Length: 3 hours                                                                                      |
| <b>Gene ID</b>              | 22925                                                                                                                                                                                                           |
| <b>Gene Symbol</b>          | PLA2R1                                                                                                                                                                                                          |
| <b>Species</b>              | Human                                                                                                                                                                                                           |
| <b>Kit Components</b>       | Pre-coated 96T strip plate, Plate sealer for 96 wells, Standard, Standard Diluent, Detection Solution A, Detection Solution B, TMB Substrate, Stop Solution, Wash Buffer (30 x concentrate), Instruction manual |
| <b>Standard Curve Range</b> | 0.156 - 10 ng/mL (example only; lot dependent)                                                                                                                                                                  |
| <b>Sensitivity</b>          | 0.065 ng/mL (example only; lot dependent)                                                                                                                                                                       |
| <b>Inter Assay</b>          | %CV < 12 (example only; lot dependent)                                                                                                                                                                          |
| <b>Intra Assay</b>          | %CV < 10 (example only; lot dependent)                                                                                                                                                                          |
| <b>Assay Type</b>           | Sandwich ELISA                                                                                                                                                                                                  |
| <b>Suitable Sample Type</b> | Serum, plasma, tissue homogenates and other biological fluids                                                                                                                                                   |
| <b>Sample Volume</b>        | 100 uL                                                                                                                                                                                                          |

**Images**

ELISA: Human PLA2R1 - Ready-To-Use ELISA Kit (Colorimetric) [NBP3-39349] - Standard Curve Reference





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP3-39349**

---

|               |                                     |
|---------------|-------------------------------------|
| NBP1-84449PEP | PLA2R1 Recombinant Protein Antigen  |
| 210-TA-005    | TNF-alpha [Unconjugated]            |
| NBP2-50248    | PLA2R1 Antibody (12-6-5) - BSA Free |
| M6000B-1      | IL-6 [HRP]                          |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. ELISA Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-39349](http://www.novusbio.com/reviews/submit/NBP3-39349)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

