

Product Datasheet

Porcine Troponin T Type 2 (cardiac) ELISA Kit (Colorimetric) NBP3-42398

Unit Size: 1 Kit

Storage of components varies. See protocol for specific instructions.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-42398

Updated 11/14/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-42398



NBP3-42398**Porcine Troponin T Type 2 (cardiac) ELISA Kit (Colorimetric)**

Product Information	
Unit Size	1 Kit
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
Storage	Storage of components varies. See protocol for specific instructions.
Product Description	
Description	Assay time: 3 hours
Gene ID	7139
Gene Symbol	TNNT2
Species	Porcine
Kit Components	Pre-coated ready to use 96-well strip plate, Plate sealer for 96 wells, Standard, Standard Diluent, Detection Reagent A, Assay Diluent A, Detection Reagent B, Assay Diluent B, TMB Substrate, Stop Solution, Wash Buffer (30 x concentrate), Instruction manual
Standard Curve Range	15.6 - 1000 pg/mL
Sensitivity	5.5 pg/mL
Inter Assay	%CV < 12
Intra Assay	%CV < 10
Assay Type	Sandwich ELISA
Spiking Recovery	83 - 94%
Suitable Sample Type	Serum, Plasma, Tissue homogenates, Cell lysates, Cell culture supernates and other biological fluids
Sample Volume	100 uL
Product Application Details	
Applications	ELISA
Recommended Dilutions	ELISA





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-42398

NBC1-28765-0.1mg	Recombinant Human Troponin T Type 2 (cardiac) His Protein
NB200-316	CUGBP1/CELF1 Antibody (3B1)
NBP3-24557-100ul	Troponin T Type 2 (cardiac) Antibody (CL13403)
NBP2-35336	Troponin I Type 3 (cardiac) Antibody (2C9) - BSA Free

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. ELISA Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-42398

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

