

# Product Datasheet

## Canine Troponin T Type 2 (cardiac) ELISA Kit (Colorimetric) NBP3-42399

Unit Size: 1 Kit

Storage of components varies. See protocol for specific instructions.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-42399](http://www.novusbio.com/NBP3-42399)

Updated 11/14/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-42399](http://www.novusbio.com/reviews/destination/NBP3-42399)



**NBP3-42399****Canine Troponin T Type 2 (cardiac) ELISA Kit (Colorimetric)**

<b>Product Information</b>	
<b>Unit Size</b>	1 Kit
<b>Concentration</b>	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
<b>Storage</b>	Storage of components varies. See protocol for specific instructions.
<b>Product Description</b>	
<b>Description</b>	Assay time: 3 hours
<b>Gene ID</b>	7139
<b>Gene Symbol</b>	TNNT2
<b>Species</b>	Canine
<b>Kit Components</b>	Pre-coated ready to use 96-well strip plate, Plate sealer for 96 wells, Standard, Standard Diluent, Detection Reagent A, Assay Diluent A, Detection Reagent B, Assay Diluent B, TMB Substrate, Stop Solution, Wash Buffer (30 x concentrate), Instruction manual
<b>Standard Curve Range</b>	31.2 - 2000 pg/mL
<b>Sensitivity</b>	12.9 pg/mL
<b>Inter Assay</b>	%CV < 12
<b>Intra Assay</b>	%CV < 10
<b>Assay Type</b>	Sandwich ELISA
<b>Spiking Recovery</b>	92 - 107%
<b>Suitable Sample Type</b>	Serum, Plasma, Tissue homogenates, Cell lysates, Cell culture supernates and other biological fluids
<b>Sample Volume</b>	100 uL
<b>Product Application Details</b>	
<b>Applications</b>	ELISA
<b>Recommended Dilutions</b>	ELISA



### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-42399**

NBC1-28765-0.1mg	Recombinant Human Troponin T Type 2 (cardiac) His Protein
NB200-316	CUGBP1/CELF1 Antibody (3B1)
NBP3-24557-100ul	Troponin T Type 2 (cardiac) Antibody (CL13403)
NBP2-35336	Troponin I Type 3 (cardiac) Antibody (2C9) - BSA Free

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. ELISA Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-42399](http://www.novusbio.com/reviews/submit/NBP3-42399)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

