

Product Datasheet

TEL2 Antibody (8B6) - BSA Free

NBP3-43043

Unit Size: 100 ug

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-43043

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-43043



NBP3-43043

TEL2 Antibody (8B6) - BSA Free

Product Information	
Unit Size	100 ug
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	8B6
Preservative	0.02% Sodium Azide
Isotype	IgG2a
Purity	Protein A or G purified
Buffer	PBS (pH 7.3), 50% glycerol
Target Molecular Weight	92 kDa


Product Description	
Description	Novus Biologicals Mouse TEL2 Antibody (8B6) - BSA Free (NBP3-43043) is a monoclonal antibody validated for use in IHC, WB, ELISA, ICC/IF and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	9894
Gene Symbol	TELO2
Species	Human
Immunogen	TEL2, telomere maintenance 2, homolog(S. cerevisiae) - (Uniprot# Q9Y4R8)

Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunoprecipitation
Recommended Dilutions	Western Blot 1:1000 - 1:4000, ELISA, Immunohistochemistry 1:20 - 1:200, Immunocytochemistry/ Immunofluorescence 1:10 - 1:100, Immunoprecipitation 1:500 - 1:1000




Images


Western Blot: TEL2 Antibody (8B6) [NBP3-43043] - A431 cells were subjected to SDS PAGE followed by western blot with (TEL2 antibody) at dilution of 1:500



Immunoprecipitation: TEL2 Antibody (8B6) [NBP3-43043] - IP Result of anti-TEL2 (IP:FNab08592, 4ug; Detection:FNab08592 1:500) with HepG2 cells lysate 2900ug.



Immunohistochemistry: TEL2 Antibody (8B6) [NBP3-43043] - Immunohistochemistry of paraffin-embedded human kidney tissue slide using (TEL2 Antibody) at dilution of 1:200





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-43043

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-96778	Mouse IgG2a Isotype Control (M2A)
NBP1-81762PEP	TEL2 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-43043

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

