

Product Datasheet

Integrin alpha 5 beta 1 Antibody (SG/19) - Chimeric - BSA Free NBP3-43579

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-43579

Updated 5/1/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-43579



NBP3-43579

Integrin alpha 5 beta 1 Antibody (SG/19) - Chimeric - BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	SG/19
Preservative	0.02% Proclin 300
Isotype	IgG
Purity	Protein A purified
Buffer	PBS
Product Description	
Description	This chimeric rabbit antibody was made using the variable domain sequences of the original Mouse IgG1 format, for improved compatibility with existing reagents, assays and techniques.
Host	Rabbit
Gene ID	3678
Gene Symbol	ITGA5
Species	Human
Specificity/Sensitivity	This antibody binds human integrin alpha 5 beta 1 near the N-terminus. Epitope mapping of SG/19 revealed that Thr82 in the 1 subunit, located at the outer face of the boundary between the I-like and hybrid domains, was the key binding determinant for this antibody. Integrin alpha5 beta 1 is also known as fibronectin receptor and binds matrix macromolecules and proteinases and thereby stimulates angiogenesis. Integrin alpha5beta1, is a member of the integrin family, is a heterodimer composed by two subunits, alpha5 (ITGA5/CD49e) and beta1 (ITGB1/CD29), and both are necessary for complete biological functions. The interaction of VLA-5 with fibronectin plays an important role in regulating inflammatory cytokine production by human articular chondrocytes (From the Cell Migration Gateway ITGA5 ITGB1
Immunogen	The original antibody was generated by immunizing BALB/c mice with Nalm-6 cells. (Uniprot# P05556)
Product Application Details	
Applications	Flow Cytometry, Functional, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation
Recommended Dilutions	Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation, Functional





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-43579

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
210-TA-005	TNF-alpha [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-43579

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

