

Product Datasheet

KRAS Antibody (D113) [DyLight 550] - Mutant NBP3-43583R

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-43583R

Updated 3/24/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-43583R



NBP3-43583R

KRAS Antibody (D113) [DyLight 550] - Mutant

| Product Information | |
|-----------------------------|---|
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | D113 |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG |
| Conjugate | DyLight 550 |
| Purity | Protein A purified |
| Buffer | 50mM Sodium Borate |
| Product Description | |
| Host | Rabbit |
| Gene ID | 3845 |
| Gene Symbol | KRAS |
| Species | Human |
| Specificity/Sensitivity | This antibody binds the G12D mutant form of human KRAS protein. This ras protein binds GDP/GTP and possess intrinsic GTPase activity, that results in conversion of nucleotide guanosine triphosphate (GTP) into guanosine diphosphate (GDP). This protein plays an important role in many signal transduction pathways and regulation of cell proliferation. Somatic KRAS mutations are highly observed in leukemias, colorectal cancer, pancreatic cancer and lung cancer. KRAS mutation is G12D which is estimated to be present in up to 37% pancreatic cancers and over 12% of colorectal cancers. |
| Immunogen | The original antibody was generated by immunizing mice with a synthetic dodecapeptide 'KLVVVGADGVGK'. (Uniprot# P01116-2) |
| Notes | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries. |
| Product Application Details | |
| Applications | ELISA |
| Recommended Dilutions | ELISA |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-43583R

| | |
|--------------------|---|
| NBP2-24891R | Rabbit IgG Isotype Control [DyLight 550] |
| H00003845-P01-10ug | Recombinant Human KRAS GST (N-Term) Protein |
| 210-TA-005 | TNF-alpha [Unconjugated] |
| NBP1-45274-0.1mg | Recombinant Human KRAS His Protein |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-43583R

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

