

Product Datasheet

PI 3-Kinase p110 beta/PIK3CB Antibody (CL9518) [PerCP] NBP3-43948PCP

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-43948PCP

Updated 3/30/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-43948PCP



NBP3-43948PCP

PI 3-Kinase p110 beta/PIK3CB Antibody (CL9518) [PerCP]

| Product Information | |
|---------------------|--|
| Unit Size | 0.1 ml |
| Concentration | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | CL9518 |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG1 |
| Conjugate | PerCP |
| Purity | Protein A purified |
| Buffer | PBS |

| Product Description | |
|---------------------|---|
| Host | Mouse |
| Gene ID | 5291 |
| Gene Symbol | PIK3CB |
| Species | Human |
| Immunogen | This antibody was generated using a synthetic peptide of P42338, with the exact immunogen sequence remaining proprietary. |

| Product Application Details | |
|-----------------------------|--|
| Applications | Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin |
| Recommended Dilutions | Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-43948PCP

| | |
|--------------------|---|
| NBP1-97005PCP | Mouse IgG1 Isotype Control (MG1) [PerCP] |
| H00005291-Q01-10ug | Recombinant Human PI 3-Kinase p110 beta/PIK3CB GST (N-Term) Protein |
| 236-EG-200 | EGF [Unconjugated] |
| NBP3-39336 | Human PI 3-Kinase p110 beta/PIK3CB ELISA Kit (Colorimetric) |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-43948PCP

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

