

Product Datasheet

hnRNP C1 + C2 Antibody (CL2593) - Azide and BSA Free NBP3-44054

Unit Size: 100 ug

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-44054

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-44054



NBP3-44054

hnRNP C1 + C2 Antibody (CL2593) - Azide and BSA Free

Product Information	
Unit Size	100 ug
Concentration	LYOPH mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	CL2593
Preservative	No Preservative
Reconstitution Instructions	Centrifuge the vial of lyophilized antibody at 12,000 x g for 20 seconds. Reconstitute by adding sterile, distilled water to achieve a final antibody concentration of 1mg/ml.
Isotype	IgG2a
Purity	Protein A purified
Buffer	Lyophilized from a 0.2 um filtered solution in PBS with Trehalose


Product Description	
Description	Novus Biologicals Mouse hnRNP C1 + C2 Antibody (CL2593) - Azide and BSA Free (NBP2-36776) is a monoclonal antibody validated for use in IHC, WB and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	3183
Gene Symbol	HNRNPC
Species	Human
Immunogen	This antibody was generated using a recombinant protein sequence of P07910, with the exact immunogen sequence remaining proprietary.

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence
Recommended Dilutions	Western Blot 1 ug/ml, Immunocytochemistry/ Immunofluorescence 2-10 ug/ml, Immunohistochemistry-Paraffin 1:1000 - 1:2500
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended. For Immunocytochemistry/ Immunofluorescence/IF, Fixation/Permeabilization with PFA/Triton X-100 is recommended.




Images


Staining in U2OS cell line with Anti-hnRNP C1 + C2 monoclonal antibody) staining in green. Microtubule- and nuclear probes are visualized in red and blue respectively (where available).




Staining in U251 cell line with Anti-hnRNP C1 + C2 monoclonal antibody) staining in green. Microtubule- and nuclear probes are visualized in red and blue respectively (where available).




Lane 1: Marker [kDa]
Lane 2: Human cell line U-251 MG




Staining in MCF7 cell line with Anti-hnRNP C1 + C2 monoclonal antibody) staining in green. Microtubule- and nuclear probes are visualized in red and blue respectively (where available).




Staining in HeLa cell line with Anti-hnRNP C1 + C2 monoclonal antibody) staining in green. Microtubule- and nuclear probes are visualized in red and blue respectively (where available).




Staining in A431 cell line with Anti-hnRNP C1 + C2 monoclonal antibody) staining in green. Microtubule- and nuclear probes are visualized in red and blue respectively (where available).




Staining of human fallopian tube shows strong nuclear immunoreactivity.




Staining of human kidney shows strong nuclear immunoreactivity.




Staining of human prostate shows strong nuclear positivity in glandular and connective tissue cells.




Staining of human rectum shows strong nuclear immunoreactivity in glandular and lamina propria cells.




Staining of human small intestine shows strong nuclear positivity in glandular epithelium and in lamina propria cells.



Staining of human testis shows nuclear positivity in seminiferous tubules.



Analysis in U-251MG cells transfected with control siRNA, target specific siRNA probe #1 and #2, using Anti-hnRNP C1 + C2 antibody. Remaining relative intensity is presented. Loading control: Anti-PPIB.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-44054

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-96778	Mouse IgG2a Isotype Control (M2A)
NBP1-51035-0.01mg	Recombinant Human hnRNP C1 + C2 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-44054

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

