

# Product Datasheet

## Cyclophilin B Antibody (CL3901) [DyLight 350] NBP3-44066UV

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-44066UV](http://www.novusbio.com/NBP3-44066UV)

Updated 3/30/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-44066UV](http://www.novusbio.com/reviews/destination/NBP3-44066UV)



**NBP3-44066UV**

Cyclophilin B Antibody (CL3901) [DyLight 350]

| Product Information |                                                                                                                  |
|---------------------|------------------------------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                                           |
| Concentration       | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                                         |
| Clonality           | Monoclonal                                                                                                       |
| Clone               | CL3901                                                                                                           |
| Preservative        | 0.05% Sodium Azide                                                                                               |
| Isotype             | IgG2a                                                                                                            |
| Conjugate           | DyLight 350                                                                                                      |
| Purity              | Protein A purified                                                                                               |
| Buffer              | 50mM Sodium Borate                                                                                               |

| Product Description |                                                                                                                           |
|---------------------|---------------------------------------------------------------------------------------------------------------------------|
| Host                | Mouse                                                                                                                     |
| Gene ID             | 5479                                                                                                                      |
| Gene Symbol         | PPIB                                                                                                                      |
| Species             | Human, Mouse, Rat                                                                                                         |
| Immunogen           | This antibody was generated using a synthetic peptide of P23284, with the exact immunogen sequence remaining proprietary. |
| Notes               | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                         |

| Product Application Details |                                                                        |
|-----------------------------|------------------------------------------------------------------------|
| Applications                | Western Blot, Immunohistochemistry-Paraffin                            |
| Recommended Dilutions       | Western Blot, Immunohistochemistry-Paraffin                            |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP3-44066UV**

---

|                 |                                             |
|-----------------|---------------------------------------------|
| NBC1-18424      | Recombinant Human Cyclophilin B Protein     |
| NBP2-59665-20ug | Recombinant Mouse Cyclophilin B His Protein |
| NB600-408-0.1mg | Collagen I Antibody - BSA Free              |
| AF972           | EMMPRIN/CD147 Antibody [Unconjugated]       |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-44066UV](http://www.novusbio.com/reviews/submit/NBP3-44066UV)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

