

# Product Datasheet

## Calbindin D-28K Antibody (CL15059) [HRP] NBP3-44160H

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-44160H](http://www.novusbio.com/NBP3-44160H)

Updated 3/30/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-44160H](http://www.novusbio.com/reviews/destination/NBP3-44160H)



**NBP3-44160H**

Calbindin D-28K Antibody (CL15059) [HRP]

**Product Information**

|                      |                                                                                                                  |
|----------------------|------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>     | 0.1 ml                                                                                                           |
| <b>Concentration</b> | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>       | Store at 4C in the dark.                                                                                         |
| <b>Clonality</b>     | Monoclonal                                                                                                       |
| <b>Clone</b>         | CL15059                                                                                                          |
| <b>Preservative</b>  | No Preservative                                                                                                  |
| <b>Isotype</b>       | IgG1                                                                                                             |
| <b>Conjugate</b>     | HRP                                                                                                              |
| <b>Purity</b>        | Protein A purified                                                                                               |
| <b>Buffer</b>        | PBS                                                                                                              |

**Product Description**

|                    |                                                                                                                                      |
|--------------------|--------------------------------------------------------------------------------------------------------------------------------------|
| <b>Host</b>        | Mouse                                                                                                                                |
| <b>Gene ID</b>     | 793                                                                                                                                  |
| <b>Gene Symbol</b> | CALB1                                                                                                                                |
| <b>Species</b>     | Human, Mouse, Rat                                                                                                                    |
| <b>Immunogen</b>   | This antibody was generated using a recombinant protein sequence of P05937, with the exact immunogen sequence remaining proprietary. |

**Product Application Details**

|                              |                                                                        |
|------------------------------|------------------------------------------------------------------------|
| <b>Applications</b>          | Western Blot, Immunohistochemistry-Paraffin                            |
| <b>Recommended Dilutions</b> | Western Blot, Immunohistochemistry-Paraffin                            |
| <b>Application Notes</b>     | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP3-44160H**

---

|            |                                           |
|------------|-------------------------------------------|
| NBP1-97074 | Mouse IgG1 Isotype Control (MG1) [HRP]    |
| NBP1-30222 | Recombinant Human Calbindin D-28K Protein |
| DBD00      | BDNF [HRP]                                |
| NBL1-08640 | Calbindin D-28K Overexpression Lysate     |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-44160H](http://www.novusbio.com/reviews/submit/NBP3-44160H)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

