

# Product Datasheet

## **Integrin alpha 4 beta 7/LPAM-1 Antibody (CL11293)** **[DyLight 350]** **NBP3-44276UV**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-44276UV](http://www.novusbio.com/NBP3-44276UV)

Updated 3/27/2025 v.20.1

**Earn rewards for product reviews and publications.**

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-44276UV](http://www.novusbio.com/reviews/destination/NBP3-44276UV)



**NBP3-44276UV**

Integrin alpha 4 beta 7/LPAM-1 Antibody (CL11293) [DyLight 350]

| Product Information |                                                                                                                  |
|---------------------|------------------------------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                                           |
| Concentration       | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                                         |
| Clonality           | Monoclonal                                                                                                       |
| Clone               | CL11293                                                                                                          |
| Preservative        | 0.05% Sodium Azide                                                                                               |
| Isotype             | IgG1                                                                                                             |
| Conjugate           | DyLight 350                                                                                                      |
| Purity              | Protein A purified                                                                                               |
| Buffer              | 50mM Sodium Borate                                                                                               |

| Product Description |                                                                                                                           |
|---------------------|---------------------------------------------------------------------------------------------------------------------------|
| Host                | Mouse                                                                                                                     |
| Gene ID             | 3676                                                                                                                      |
| Gene Symbol         | ITGA4, ITGA4, Itga5, Itga5                                                                                                |
| Species             | Human                                                                                                                     |
| Immunogen           | This antibody was generated using a synthetic peptide of P13612, with the exact immunogen sequence remaining proprietary. |
| Notes               | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                         |

| Product Application Details |                                                                        |
|-----------------------------|------------------------------------------------------------------------|
| Applications                | Western Blot                                                           |
| Recommended Dilutions       | Western Blot                                                           |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-44276UV**

---

|              |                                                |
|--------------|------------------------------------------------|
| NBP1-97005UV | Mouse IgG1 Isotype Control (MG1) [DyLight 350] |
| 210-TA-005   | TNF-alpha [Unconjugated]                       |
| 5397-A3-050  | Integrin alpha 4 beta 7/LPAM-1 [Unconjugated]  |
| AF3628       | CD31/PECAM-1 Antibody [Unconjugated]           |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-44276UV](http://www.novusbio.com/reviews/submit/NBP3-44276UV)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

