

# Product Datasheet

## CD19 Antibody (CL15454) [DyLight 550] NBP3-44277R

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-44277R](http://www.novusbio.com/NBP3-44277R)

Updated 3/27/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-44277R](http://www.novusbio.com/reviews/destination/NBP3-44277R)



**NBP3-44277R**

CD19 Antibody (CL15454) [DyLight 550]

| <b>Product Information</b> |  |
|----------------------------|--|
| <b>Unit Size</b>           | 0.1 ml   |
| <b>Concentration</b>       | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>             | Store at 4C in the dark.   |
| <b>Clonality</b>           | Monoclonal   |
| <b>Clone</b>               | CL15454  |
| <b>Preservative</b>        | 0.05% Sodium Azide   |
| <b>Isotype</b>             | IgG2b  |
| <b>Conjugate</b>           | DyLight 550  |
| <b>Purity</b>              | Protein A purified   |
| <b>Buffer</b>              | 50mM Sodium Borate   |

| <b>Product Description</b> |   |
|----------------------------|---|
| <b>Host</b>                | Mouse   |
| <b>Gene ID</b>             | 930   |
| <b>Gene Symbol</b>         | CD19  |
| <b>Species</b>             | Human   |
| <b>Immunogen</b>           | This antibody was generated using a synthetic peptide of P15391, with the exact immunogen sequence remaining proprietary. |
| <b>Notes</b>               | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.   |

| <b>Product Application Details</b> |  |
|------------------------------------|--|
| <b>Applications</b>                | Western Blot, Immunohistochemistry-Paraffin                            |
| <b>Recommended Dilutions</b>       | Western Blot, Immunohistochemistry-Paraffin                            |
| <b>Application Notes</b>           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-44277R**

---

|                    |  |
|--------------------|--|
| NBP2-27231R        | Mouse IgG2b Isotype Control (MPC-11) [DyLight 550] |
| H00000930-P01-10ug | Recombinant Human CD19 GST (N-Term) Protein        |
| 210-TA-005         | TNF-alpha [Unconjugated]                           |
| 9269-CD-050        | CD19 [Unconjugated]                                |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-44277R](http://www.novusbio.com/reviews/submit/NBP3-44277R)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

