

# Product Datasheet

## SSEA-1 Antibody (CL6334) [Biotin] NBP3-44428B

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-44428B](http://www.novusbio.com/NBP3-44428B)

Updated 3/30/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-44428B](http://www.novusbio.com/reviews/destination/NBP3-44428B)



**NBP3-44428B**

SSEA-1 Antibody (CL6334) [Biotin]

| <b>Product Information</b> |                                                                                                                  |
|----------------------------|------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>           | 0.1 ml                                                                                                           |
| <b>Concentration</b>       | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>             | Store at 4C in the dark.                                                                                         |
| <b>Clonality</b>           | Monoclonal                                                                                                       |
| <b>Clone</b>               | CL6334                                                                                                           |
| <b>Preservative</b>        | 0.05% Sodium Azide                                                                                               |
| <b>Isotype</b>             | IgG2a                                                                                                            |
| <b>Conjugate</b>           | Biotin                                                                                                           |
| <b>Purity</b>              | Protein A purified                                                                                               |
| <b>Buffer</b>              | PBS                                                                                                              |

| <b>Product Description</b> |                                                                                                                           |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------|
| <b>Host</b>                | Mouse                                                                                                                     |
| <b>Gene ID</b>             | 2526                                                                                                                      |
| <b>Gene Symbol</b>         | FUT4                                                                                                                      |
| <b>Species</b>             | Human                                                                                                                     |
| <b>Immunogen</b>           | This antibody was generated using a synthetic peptide of P22083, with the exact immunogen sequence remaining proprietary. |

| <b>Product Application Details</b> |                                                                        |
|------------------------------------|------------------------------------------------------------------------|
| <b>Applications</b>                | Western Blot                                                           |
| <b>Recommended Dilutions</b>       | Western Blot                                                           |
| <b>Application Notes</b>           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-44428B**

---

|                    |                                               |
|--------------------|-----------------------------------------------|
| IC003B             | Mouse IgG2a Isotype Control (20102) [Biotin]  |
| H00002526-Q01-10ug | Recombinant Human SSEA-1 GST (N-Term) Protein |
| DC140              | CD14 [HRP]                                    |
| NBP2-10793         | SSEA-1 Overexpression Lysate                  |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-44428B](http://www.novusbio.com/reviews/submit/NBP3-44428B)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

