

# Product Datasheet

## UBE3B Antibody (OTI3F10) [DyLight 755] NBP3-44687IR

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-44687IR](http://www.novusbio.com/NBP3-44687IR)

Updated 5/20/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-44687IR](http://www.novusbio.com/reviews/destination/NBP3-44687IR)



**NBP3-44687IR**

UBE3B Antibody (OTI3F10) [DyLight 755]

| Product Information         |                                                                                                                         |
|-----------------------------|-------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                  |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                             |
| Storage                     | Store at 4C in the dark.                                                                                                |
| Clonality                   | Monoclonal                                                                                                              |
| Clone                       | OTI3F10                                                                                                                 |
| Preservative                | 0.05% Sodium Azide                                                                                                      |
| Isotype                     | IgG2a                                                                                                                   |
| Conjugate                   | DyLight 755                                                                                                             |
| Purity                      | >95%                                                                                                                    |
| Buffer                      | 50mM Sodium Borate                                                                                                      |
| Product Description         |                                                                                                                         |
| Host                        | Mouse                                                                                                                   |
| Gene ID                     | 89910                                                                                                                   |
| Gene Symbol                 | UBE3B                                                                                                                   |
| Species                     | Human, Mouse, Rat                                                                                                       |
| Immunogen                   | Human recombinant protein fragment corresponding to amino acids 768-1068 of human UBE3B (NP_569733) produced in E.coli. |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                       |
| Product Application Details |                                                                                                                         |
| Applications                | Immunohistochemistry                                                                                                    |
| Recommended Dilutions       | Immunohistochemistry                                                                                                    |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                  |



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to NBP3-44687IR**

---

|                |                                          |
|----------------|------------------------------------------|
| NBP1-92559PEP  | UBE3B Recombinant Protein Antigen        |
| NB300-581-10uL | SERCA2 ATPase Antibody (2A7-A1)          |
| NBP2-22106     | PPAR gamma/NR1C3 Antibody - BSA Free     |
| H00004595-M01  | MYH Antibody (4D10) - Azide and BSA Free |

---

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-44687IR](http://www.novusbio.com/reviews/submit/NBP3-44687IR)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

