

# Product Datasheet

## BAP1 Antibody (BAP1/13501R) NBP3-45077

Unit Size: 100 ug

Store at -20 to -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-45077](http://www.novusbio.com/NBP3-45077)

Updated 3/19/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-45077](http://www.novusbio.com/reviews/destination/NBP3-45077)



**NBP3-45077**

BAP1 Antibody (BAP1/13501R)

| Product Information     |                                                 |
|-------------------------|-------------------------------------------------|
| Unit Size               | 100 ug                                          |
| Concentration           | 1 mg/ml                                         |
| Storage                 | Store at -20 to -80C. Avoid freeze-thaw cycles. |
| Clonality               | Monoclonal                                      |
| Clone                   | BAP1/13501R                                     |
| Preservative            | No Preservative                                 |
| Isotype                 | IgG Kappa                                       |
| Purity                  | Protein A or G purified                         |
| Buffer                  | 10 mM PBS                                       |
| Target Molecular Weight | 91 kDa                                          |

| Product Description |                                                                                                                                                                                                                                                                                                  |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description         | 1.0 mg/ml of Ab produced in CHO cell mammalian-based expression system.. Prepared in 10 mM PBS WITHOUT BSA & azide. Also available at 200 ug/ml with BSA & azide.<br><br>Antibody with azide - store at 2 to 8C. Antibody without azide - store at -20 to -80C. Non-hazardous. No MSDS required. |
| Host                | Rabbit                                                                                                                                                                                                                                                                                           |
| Gene ID             | 8314                                                                                                                                                                                                                                                                                             |
| Gene Symbol         | BAP1                                                                                                                                                                                                                                                                                             |
| Species             | Human                                                                                                                                                                                                                                                                                            |
| Immunogen           | Recombinant fragment (around aa400-729) of human BAP1 protein (Uniprot: Q92560). Exact sequence is proprietary.                                                                                                                                                                                  |

| Product Application Details |                                                                                                                                                                                                         |
|-----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Applications                | Western Blot, Immunohistochemistry                                                                                                                                                                      |
| Recommended Dilutions       | Western Blot : 2-4ug/ml, Immunohistochemistry : 1-2ug/ml                                                                                                                                                |
| Application Notes           | Immunohistochemistry: 30 min at RT. Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95C followed by cooling at RT for 20 minutes. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP3-45077**

---

|               |                                                     |
|---------------|-----------------------------------------------------|
| HAF008        | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160        | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-55462PEP | BAP1 Recombinant Protein Antigen                    |
| NBL1-07915    | BAP1 Overexpression Lysate                          |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-45077](http://www.novusbio.com/reviews/submit/NBP3-45077)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

