

Product Datasheet

TRAPPC3 Antibody - BSA Free

NBP3-45740

Unit Size: 100 ug

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-45740

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-45740



NBP3-45740

TRAPPC3 Antibody - BSA Free

Product Information

Unit Size	100 ug
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS, 50% Glycerol, pH (7.3)

Product Description

Description	Novus Biologicals Rabbit TRAPPC3 Antibody - BSA Free (NBP3-45740) is a polyclonal antibody validated for use in IHC, WB, ELISA, ICC/IF and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	27095
Gene Symbol	TRAPPC3
Species	Human, Mouse
Immunogen	trafficking protein particle complex 3. (Uniprot# O43617)


Product Application Details

Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunoprecipitation
Recommended Dilutions	Western Blot 1:200-1:2000, ELISA, Immunohistochemistry 1:20-1:200, Immunocytochemistry/ Immunofluorescence 1:10-1:100, Immunoprecipitation 1:200-1:2000
Application Notes	Theoretical Molecular Weight is 20 kDa, 30 kDa




Images


Immunoprecipitation: TRAPPC3 Antibody - BSA Free [NBP3-45740] - IP Result of anti-TRAPPC3 (IP: 3ug; Detection: 1:800) with mouse liver tissue lysate 4000ug.



Immunohistochemistry: TRAPPC3 Antibody - BSA Free [NBP3-45740] - Immunohistochemistry of paraffin-embedded human placenta using (TRAPPC3 antibody) at dilution of 1:100



Western Blot: TRAPPC3 Antibody - BSA Free [NBP3-45740] - Mouse liver tissue were subjected to SDS PAGE followed by western blot with (TRAPPC3 antibody) at dilution of 1:800





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-45740

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP1-99070-100ug	Recombinant Human TRAPPC3 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-45740

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

