

Product Datasheet

PHLPP Antibody - BSA Free

NBP3-47157

Unit Size: 100 ug

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-47157

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-47157



NBP3-47157


PHLPP Antibody - BSA Free

Product Information	
Unit Size	100 ug
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS, 50% Glycerol, pH (7.3)
Target Molecular Weight	185 kDa
Product Description	
Description	Novus Biologicals Rabbit PHLPP Antibody - BSA Free (NBP3-47157) is a polyclonal antibody validated for use in IHC, WB, ELISA, ICC/IF and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	23239
Gene Symbol	PHLPP1
Species	Human, Mouse
Immunogen	PH domain and leucine rich repeat protein phosphatase. (Uniprot# O60346)
Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunoprecipitation
Recommended Dilutions	Western Blot 1:500-1:2000, ELISA, Immunohistochemistry 1:20-1:200, Immunocytochemistry/ Immunofluorescence 1:20-1:200, Immunoprecipitation 1:500-1:2000



Images


Immunohistochemistry: PHLPP Antibody - BSA Free [NBP3-47157] - Immunohistochemistry of paraffin-embedded human liver slide using (PHLPP Antibody) at dilution of 1:50




Immunocytochemistry/ Immunofluorescence: PHLPP Antibody - BSA Free [NBP3-47157] - Immunofluorescent analysis of HepG2 cells, using PHLPP antibody at 1:50 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).



Immunoprecipitation: PHLPP Antibody - BSA Free [NBP3-47157] - IP Result of anti-PHLPP (IP: , 4ug; Detection: 1:1000) with mouse brain tissue lysate 5000ug.



Western Blot: PHLPP Antibody - BSA Free [NBP3-47157] - Mouse brain tissue were subjected to SDS PAGE followed by western blot with (PHLPP antibody) at dilution of 1:1000





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-47157

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-33887PEP	PHLPP Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-47157

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

