

Product Datasheet

DNAJB4 Antibody - BSA Free

NBP3-47990

Unit Size: 100 ug

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-47990

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-47990



NBP3-47990

DNAJB4 Antibody - BSA Free

Product Information

Unit Size	100 ug
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS, 50% Glycerol, pH (7.3)
Target Molecular Weight	38 kDa

Product Description

Description	Novus Biologicals Rabbit DNAJB4 Antibody - BSA Free (NBP3-47990) is a polyclonal antibody validated for use in IHC, WB, ELISA and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	11080
Gene Symbol	DNAJB4
Species	Human, Mouse, Rat
Immunogen	DnaJ(Hsp40) homolog, subfamily B, member 4. (Uniprot# Q9UDY4)


Product Application Details

Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1:1000-1:4000, ELISA, Immunohistochemistry 1:20-1:200, Immunocytochemistry/ Immunofluorescence 1:20-1:200




Images


Immunohistochemistry: DNAJB4 Antibody - BSA Free [NBP3-47990] - Immunohistochemistry of paraffin-embedded human pancreas cancer using (DNAJB4 antibody) at dilution of 1:100



Western Blot: DNAJB4 Antibody - BSA Free [NBP3-47990] - HeLa cells were subjected to SDS PAGE followed by western blot with (DNAJB4 Antibody) at dilution of 1:2000



Immunocytochemistry/ Immunofluorescence: DNAJB4 Antibody - BSA Free [NBP3-47990] - Immunofluorescent analysis of MCF-7 cells, using (DNAJB4 antibody) at 1:50 dilution and Rhodamine-labeled goat anti-rabbit IgG (red)





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-47990

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-23005	Recombinant Human DNAJB4 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-47990

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

