

Product Datasheet

HAPLN2 Antibody (GT107)

NBP3-48853

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-48853

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-48853



NBP3-48853

HAPLN2 Antibody (GT107)

Product Information	
Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	GT107
Preservative	No Preservative
Isotype	IgG2b
Purity	Protein A purified
Buffer	PBS


Product Description	
Description	Novus Biologicals Mouse HAPLN2 Antibody (GT107) (NBP3-48853) is a recombinant monoclonal antibody validated for use in IHC and WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	60484
Gene Symbol	HAPLN2
Species	Mouse
Immunogen	The immunogen used to generate this antibody corresponds to human HAPLN2.

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry
Recommended Dilutions	Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin




Images


Various tissue extracts (50 ug) were separated by 10% SDS-PAGE, and the membrane was blotted with HAPLN2 antibody [GT107] (NBP3-48853) diluted at 1:1000. The HRP-conjugated anti-mouse IgG antibody was used to detect the primary antibody. Corresponding RNA expression data are based on NCBI database.



HAPLN2 antibody [GT107] detects HAPLN2 protein by immunohistochemical analysis. Sample: Paraffin-embedded mouse striatum. HAPLN2 stained by HAPLN2 antibody [GT107] (NBP3-48853) diluted at 1:200. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



HAPLN2 antibody [GT107] detects HAPLN2 protein by immunohistochemical analysis. Sample: Paraffin-embedded mouse tissues. HAPLN2 stained by HAPLN2 antibody [GT107] (NBP3-48853) diluted at 1:200. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-48853

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP2-27231	Mouse IgG2b Isotype Control (MPC-11)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-48853

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

