

Product Datasheet

EphA3 Antibody (IIIA4) [DyLight 755] - Chimeric NBP3-48887IR

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-48887IR

Updated 9/2/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-48887IR



NBP3-48887IR

EphA3 Antibody (IIIA4) [DyLight 755] - Chimeric

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	IIIA4
Preservative	0.05% Sodium Azide
Isotype	IgG Kappa
Conjugate	DyLight 755
Purity	Protein A purified
Buffer	50mM Sodium Borate

Product Description	
Host	Rabbit
Gene ID	2042
Species	Human, Mouse
Specificity/Sensitivity	The antibody binds to the native EphA3 globular ephrin-binding domain. The antibody targets a site closely adjacent to the heterotetramerization site on the N-terminal of EPHA3's extracellular domain adjacent to the ligand-binding site and has a high affinity for EPHA3.
Immunogen	The original antibody was prepared by intraperitoneal immunization of a BALB/c mouse with LK63 cells (Uniprot#: P29320)
Notes	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.

Product Application Details	
Applications	Western Blot, ELISA, Flow Cytometry, Immunocytochemistry/Immunofluorescence, In vivo assay, Immunoprecipitation, Surface Plasmon Resonance
Recommended Dilutions	Western Blot, Flow Cytometry, ELISA, Immunocytochemistry/Immunofluorescence, Immunoprecipitation, In vivo assay, Surface Plasmon Resonance
Application Notes	Optimal dilution of this antibody should be experimentally determined.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-48887IR

NBP3-17708PEP	EphA3 Recombinant Protein Antigen
236-EG-200	EGF [Unconjugated]
640-A3-200	EphA3
210-TA-005	TNF-alpha [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-48887IR

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

