

# Product Datasheet

## HIST1H2AK Antibody (RB26883) [DyLight 405] NBP3-52282V

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-52282V](http://www.novusbio.com/NBP3-52282V)

Updated 7/20/2026 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-52282V](http://www.novusbio.com/reviews/destination/NBP3-52282V)



**NBP3-52282V**

HIST1H2AK Antibody (RB26883) [DyLight 405]

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                    |
| Clonality           | Polyclonal                                                                                  |
| Clone               | RB26883                                                                                     |
| Preservative        | 0.05% Sodium Azide                                                                          |
| Isotype             | IgG                                                                                         |
| Conjugate           | DyLight 405                                                                                 |
| Purity              | Antigen and protein A Affinity-purified                                                     |
| Buffer              | 50mM Sodium Borate                                                                          |

| Product Description |                                                                                                                                                                                                                                                         |
|---------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Host                | Rabbit                                                                                                                                                                                                                                                  |
| Gene ID             | 8330                                                                                                                                                                                                                                                    |
| Gene Symbol         | H2AC15                                                                                                                                                                                                                                                  |
| Species             | Human                                                                                                                                                                                                                                                   |
| Immunogen           | This HIST1H2AK antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 1-30 amino acids from the N-terminal region of human HIST1H2AK.<br>(Accession#:NP_003500.1;NP_003501.1;NP_003502.1;NP_003505.1;NP_066408.1) |
| Notes               | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                                                                                       |

| Product Application Details |                                                                        |
|-----------------------------|------------------------------------------------------------------------|
| Applications                | Western Blot                                                           |
| Recommended Dilutions       | Western Blot                                                           |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP3-52282V**

---

|                    |                                                  |
|--------------------|--------------------------------------------------|
| NBP2-24891V        | Rabbit IgG Isotype Control [DyLight 405]         |
| H00008330-P01-10ug | Recombinant Human HIST1H2AK GST (N-Term) Protein |
| NBL1-11570         | H4/e Overexpression Lysate                       |
| NB100-56346        | Histone H2a Antibody - BSA Free                  |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-52282V](http://www.novusbio.com/reviews/submit/NBP3-52282V)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

